

minispec Sales Information 11/09: Unified minispec Plus Version V3.1



Bruker Optik

Dear all:

We are happy to announce release of the unified minispec Software with minispec Plus (m+) V3.1. This release contains new features and new applications and is compatible for all three minispec product lines, mq, mq-one and LF series. Key features are:

- All-in-One Installation¹ with minispec Plus, minispec.exe, ExpSpel etc.
- New and unique m+ GUI with PilotScan, TipTexts and Calibration Validation
- New m+ Applications: Fat in Liquor, Indirect SFC, Solid Echo, chemometrics for mq-one
- Support of OPUS 6.5 and animal weighing mode as well as multi-quant, multi-evaluation models (for BCA and chemometrical applications)
- Full Temperature Configuration with setting, profiles and ramps (together with mq-BVT)
- m+ display text in: English, German, French, Spanish, Italian, Turkish, Chinese, Japanese, Russian (can be chosen for each individual user)
- Automatic back-end upgrades
- Integrated Update Settings, integrated toggle mode of optical sample detector
- Improved data management (edit reference values, hide expired calibration)

minispec m+ V3.1 is released for Windows® XP SP2 or SP3 platform English or German version. Office language packs can be installed to offer native language support for Office (read release and installation notes for more details). In case the PC is not provided by us care has to be taken, that Microsoft® Office 2007 or Excel 2007 English version is installed on the PC prior to installing m+ V3.1.

Demo Versions are available. Please contact your regional minispec responsible.

Yours Sincerely
minispec Team

¹ Depending on the licenses purchased, the respective software packages are installed accordingly.